

Downloading Your Animas Pump


Step 1. Suspend your pump.

Step 2. Detach the pump from your body

Step 3. Place the pump in the blue stand



Step 4. Place the stand next to the transmitter (make sure the back of the pump is facing the transmitter).



Step 5. Make sure the pump screen is lit up by pressing the **Contrast Button** on the top of the pump



Step 6. The blue light on the transmitter should start flashing rapidly when the pump is downloading

Step 7. When the green light comes on, the pump is finished downloading and you can re-attach your pump and restart the insulin delivery