

Downloading your Medtronic Pump

Step 1. Go to internet explorer and go to favorites and click on Medtronic Carelink Personal

<https://carelink.minimed.eu>

The screenshot shows the CareLink Personal website homepage. At the top left is the Medtronic logo. To the right is the CareLink Personal logo with the tagline 'THERAPY MANAGEMENT SOFTWARE FOR DIABETES'. Below the logos is a welcome message: 'Welcome to the CareLink® Personal Therapy Management Software for Diabetes. This Web-based system is designed to help you take information from all of your diabetes management tools – your insulin pump, continuous glucose monitor, blood glucose meter(s), and logbook – and organize it into easy-to-read charts, graphs and tables. These reports can help you and your healthcare provider discover trends and other information that can lead to improved therapy management for greater control.' Below this is a diagram showing the flow of data from 'Medtronic Insulin Pump or Continuous Glucose Monitor', 'Blood Glucose Meter', and 'Logbook' to 'Medtronic CareLink* Software'. A box titled 'Features of the CareLink Personal Software:' lists 'Personal treatment reports with the information you need'. On the right side, there is a 'Sign Up Now' button, a 'Sign In' button, and a 'Forgot your password?' link. Below the sign in section are 'Other Resources' including 'MiniMed.com' and 'Medtronic MiniMed Online Store'.

Step 2. Sign in with your username and password

Step 3. Click on Upload tab

The screenshot shows the CareLink Personal software interface. At the top, there are four tabs: 'Home', 'Upload', 'Logbook', and 'Reports'. The 'Upload' tab is selected. Below the tabs, there is a table with three rows and two columns. The first column is empty, and the second column is labeled 'Serial #'. The three rows all contain the number '410455'. To the right of the table, there is a section titled 'What Can I Do Next?' with three options: 'Upload Data from My Device' (circled in red), 'Enter Data into My Logbook', and 'Generate Reports'.

	Serial #
	410455
	410455
	410455

Step 4. Follow prompts

Step 5. Make sure the pump remains near the USB once plugged in.



Step 6. The pump will be automatically suspended as it uploads

Step 7. Once upload is completed click on the reports tab

Step 8. Click on each of the following reports below for the last 2 weeks and print.

Quick view summary

Logbook diary

Modal day hourly

Device settings

The screenshot shows the CareLink software interface. On the left is a sidebar menu with the following items: Quick View Summary, Daily Summary, Logbook Diary, Data Table, Modal Day Periods, Modal Day Hourly, Trends Summary, Device Settings, Sensor Daily Overlay, Sensor Weekly Logbook, Sensor Overlay by Meal, and Data Export (CSV). The 'Quick View Summary' item is selected and highlighted. The main content area shows the 'Quick View Summary' report for 'Two weeks ending: 16/03/2015'. The report includes a graph of glucose and insulin levels and a table of logbook data. A description of the report is provided: 'This report shows graphical summaries of glucose and insulin along with statistical information and logbook data in tables for a two-week period. It is designed to assist your healthcare provider with a one-page summary of the most important information about your therapy.'