

24 HOUR URINE COLLECTION GUIDE

TEST	BOTTLE	CODE
Aldosterone	Plain	24ADO
Arsenic	Metal Free (Acid Washed)	24ARS
Chloride	Plain	24CL
Calcium	Plain	24CA
Cadmium	Plain	24CD
Cystine	Plain	CYSMT
Cortisol (Urine Free Cortisol)	Plain	24UFC
Creatinine	Plain	24R
Creatinine Clearance	Plain	CRCL
Catecholamines	25mL 6N HCl	24CTF
Citrate	Plain	24CIT
Dopamine	Plain	24DOP
Glucose	Plain	24GLU
Hydroxyproline	25mL of 6HCl	24HP
Immunofixation/ImmunoElectrophoresis/Bence Jones Protein	Plain	24UIF
Lytes (Na, K, Cl)	Plain	24UL
Litho Samples	Plain	“ULIT1” & “ULIT2” (enter both)
Lead	Metal Free (Acid Washed)	24PB
Magnesium	Plain	24MG
Mercury	Metal Free (Acid Washed)	24HG24HG
Free Metanephrines	25mL of 6HCl	24MET
Microalbumin/Albumin	Plain	24MAL
Oxalate	Plain	24OX
Potassium	Plain	24K
Phosphate/ Phosphorus	Plain	24PO4
Protein	Plain	24TP
Protein Electrophoresis	Plain	24UPE
Porphobilinogen	Plain (STAT – Transport Keep Cold!!)	24PBG
Porphorins	5g Sodium Carbonate	24POR
Sodium	Plain	24NA
Uric Acid	Plain	24URI
Urea Clearance	Plain	URCL
Urea	Plain	24URE
VMA	25mL of 6HCl	24VMA
Zinc	Plain	24ZN
5HIAA (5-Hydroxyindoleacetic acid)	25mL of 6HCl	24HIA

For All Send Out Urines:

1. Make sure bottles are labeled before giving it out to patients.
2. Make sure to print & attach **patient instructions** with the urine requisition
3. Enter the **height & Weight** of the patient esp with the CrClearance Urine Collection
4. Write the **start and end time** of the urine 24 hour collection
5. If the sample has to be protected from light, wrap in foil. The sample must be labeled on the bottle and also on the tin foil(after wrapping)
6. If its PG D2 & Methylhistamine urine testing, pls make sure to ask the patient to deliver it straight to CoreLAB Receiving (open 24/7) kept cold while on transport!