

# DEPARTMENT OF LABORATORY MEDICINE DIVISION OF MICROBIOLOGY

Document Name:

Collection Instructions of Stool for C.difficile Testing

Document #: 133669 Status: Current

Uncontrolled When Printed

# COLLECTION OF STOOL FOR C.difficile TESTING

### **Materials supplied:**

Sterile specimen container Biohazard transport bag Plastic pipette Requisition completed by physician or clinic

## **Materials needed but not supplied:**

A clean dry container

#### **Collection guidelines:**

Collect 1 loose watery or soft stool specimen

#### **Rejection criteria:**

Testing on formed or hard stool is **NOT** performed

### **Collection procedure:**

- 1. Pass loose watery or soft stool directly into a clean dry container. Soft stool is defined as stool assuming the shape of its container.
- 2. Open the lid of the specimen container.
- 3. Using the plastic pipette, transfer 5 10 mL ( $\sim 1 2$  teaspoon in volume) of specimen to sterile specimen container.
- 4. Discard pipette into garbage.
- 5. Screw cap tightly on the specimen container.
- 6. Wash your hands thoroughly.
- 7. Write the patient's full name, date and time of collection on the specimen container label.
- 8. Place collected specimen inside the biohazard transport bag.
- 9. Place the completed requisition into pouch on the side of the biohazard bag.
- 10. Collected specimen must be stored in the refrigerator.

Deliver specimens to the Microbiology laboratory within 24 hours of collection. The microbiology laboratory will receive specimens 24 hours/7 days a week.

Microbiology Laboratory: St. Michael's Hospital, Cardinal Carter wing, 2<sup>nd</sup> floor – Room 2044

Telephone #: 416-864-5381

Authority for Issue: Dr. M. Tadros	Authorized Date: 3/6/2019
Version: 1.0.1	Effective Date: 3/6/2019