

 St. Michael's	<p align="center">DEPARTMENT OF LABORATORY MEDICINE DIVISION OF MICROBIOLOGY</p>	<p align="right">Document Name: Collection Instructions of Stool for Ova and Parasites Document #: 133671 Status: Current</p>
--	---	---

Uncontrolled When Printed

COLLECTION OF STOOL FOR OVA AND PARASITES

Materials supplied:

Sterile transport container with SAF fixative for Ova & Parasites x3
 Biohazard transport bag x3
 Wooden Spatula x3
 Requisition completed by physician or clinic with Information on symptoms, travel, etc. (use a separate sheet of paper if necessary).

Materials needed but not supplied:

Clean wide container

Collection guidelines:

****Collect 1 specimen/day. Addition specimens should be 2 days apart****

**** A maximum of 3 specimens should be collected****

****Do not take any medicines (laxatives, anti-diarrhea drugs) containing mineral oil, barium, bismuth, kaolin, magnesium, antibiotics, etc., for at least 5 – 7 days prior to collection****

Collection Procedure:

1. Pass stool into a wide-mouthed container or onto clean paper.
NOTE: DO NOT MIX STOOL WITH URINE OR WATER FROM TOILET BOWL.
2. Using the spoon attached to the lid of the container or the wooden spatula provided, transfer a quantity of stool specimen into the container until the liquid reaches the **“FILL LINE”** on the outside of the container.
IMPORTANT: Mix the specimen well with the fluid in the container.
3. Replace the lid on the container and tighten securely.
4. Shake until mixed well.
5. Wash and dry hands thoroughly.
6. Write patient's full name, date and time the specimen was collected on container label.
7. Place specimen in the biohazard transport bag (one specimen/bag).
8. Place the completed requisition in side pouch of biohazard transport bag.
9. Keep collected specimen at room temperature.

After collection, all 3 specimens can be delivered to the Microbiology Laboratory at the same time.

Microbiology Laboratory:	St. Michael's Hospital, Cardinal Carter wing, 2nd floor – Room 2044	
Receiving Hours:	Weekdays	08:30 – 22:00
	Weekends/Holidays	08:00 – 16:00
Telephone #:	416-864-5381	

Authority for Issue: Dr. L. Matukas	Authorized Date: 7/17/2020
Version: 1.0.1	Effective Date: 7/17/2020
Any document appearing in paper form is uncontrolled and should be checked against the master electronic current version prior to use. Only original printed material with the “CONTROLLED” water mark may exist in designated locations. The controlled printed document should only be used when the electronic version is unavailable. Unauthorized photocopies or alterations of this document are uncontrolled documents.	

 St. Michael's	<p style="text-align: center;">DEPARTMENT OF LABORATORY MEDICINE DIVISION OF MICROBIOLOGY</p>	<p>Document Name: Collection Instructions of Stool for Ova and Parasites</p> <p>Document #: 133671 Status: Current</p>
--	--	--

Uncontrolled When Printed

Authority for Issue: Dr. L. Matukas	Authorized Date: 7/17/2020
Version: 1.0.1	Effective Date: 7/17/2020
<p>Any document appearing in paper form is uncontrolled and should be checked against the master electronic current version prior to use. Only original printed material with the "CONTROLLED" water mark may exist in designated locations. The controlled printed document should only be used when the electronic version is unavailable. Unauthorized photocopies or alterations of this document are uncontrolled documents.</p>	